


## **Detector Components**

## MCT SERIES CRYOGENIC PHOTODETECTOR/AMPLIFIER



## Application Note

Part No:

MCT10-0025-E-LN6N

This unit is a high performance cryogenically operated HgCdTe photodetector/amplifier. The should be unit at LN2 temperature before turning on power to the amplifier. Α funnel is provided to assist in the filling of the dewar, which best accomplished is by gradually filling and topping off over a several minute period.

The amplifier has a dual gain function controlled by a switch on the backplate. The HI (up) position is x10 above the LO (down) position. Output is thru a BNC-type cable, and power is connected thru a shielded multi-wire cable terminated in a 9-pin Dsub connector or solder leads.

SPECIFICATIONS	
Active Area	.25 mm x .25 mm
Spectral Range	2 – 12+ um
Detectivity (D*pk)	>4.0 x 10 <sup>10</sup> cm-Hz1/2/W, min
Dewar Hold Time	12 hours minimum with liquid $N_2$
Field of View	60° nominal
Responsivity (pk), at amplifier out	4.0 x 10 <sup>6</sup> V/W HI ; 10 <sup>5</sup> V/W LO
Bandwidth	5 Hz - 50kHz + typ
Connections	BNC signal coaxial cable with 3 lead shielded power cable. Red = +V, Black = -V, White/Shield = ground Note: A DB9 connector is provided on units purchased with optional PS-1 Low Noise Power Supply