


## MCT SERIES CRYOGENIC PHOTODETECTOR/AMPLIFIER



Part No: MCT14-020-E-LN6N s/n 121119

# Application Note

This unit is a high performance cryogenically operated HgCdTe photodetector/amplifier. The unit should be at LN2 temperature before turning on power to the amplifier. A funnel is provided to assist in the filling of the dewar, which is best accomplished by gradually filling and topping off over a several minute period.

The amplifier has a dual gain function controlled by a switch on the backplate. The HI (up) position is x10 above the LO (down) position. Output is thru a BNC-type cable, and power is connected thru a shielded multi-wire cable terminated in a 9-pin Dsub connector or solder leads.

SPECIFICATIONS	
Active Area	2 mm x 2 mm
Spectral Range	2 – 14+ um
Dewar window	ARGE
Dewar Hold Time	12 hours minimum with liquid N <sub>2</sub>
Field of View	60° nominal
Responsivity (pk), at amplifier out	9.6 x 10 <sup>5</sup> HI ; 10 <sup>4</sup> LO V/W
<b>Detector Resistance</b>	47 ohms
Detector Bias (set internally)	20 mA



### **Detector Components**

## MCT SERIES CRYOGENIC PHOTODETECTOR/AMPLIFIER

Noise Voltage @ 10kHz, amp out	$4.3 \times 10^{-6} \text{ HI}$ ; x $10^{-7} \text{ LO V/Hz}^{1/2}$
Detectivity (pk,10kHz,1Hz)	4.4 x 10 <sup>10</sup> cm-Hz <sup>1/2</sup> /W
Bandwidth	5 Hz - 50kHz + typ
Connections	BNC signal coaxial cable with 3 lead shielded power cable.  Red = +V, Black = -V, White/Shield = ground  Note: A DB9 connector is provided on units purchased  with optional PS-1 Low Noise Power Supply

#### Part No: LN-series Receivers

### **DB-9 PIN OUT**

- L NO CONNECT 6 +V
- 2 NO CONNECT 7 -V 3 NO CONNECT 8 GND
- NO CONNECT 9 CASE GND
- 5 NO CONNECT

Electro - Optical Systems Inc, 1288 Valley Forge Rd, Ste 49, Phoenixville, PA 19460 Tel: 610-935-5838 (Fax: -8548) www.eosystems.com sales@eosystems.com